

WHEELCHAIRS | USER GUIDE

rainbow care

BOBBY ASTRO 12 LIGHTWEIGHT WHEELCHAIR



For illustration purposes only

VERMEIREN

BOBBY ASTRO 12 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW



Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with detachable foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aircraft grade aluminium lightweight frame
- Designed in Belgium

For illustration purposes only

VERMEIREN LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



FLIP UP ARM REST

Press button to release lock and flip up arm rest



DETACHABLE FOOTREST

Pull the lever to release lock and flip up arm rest



STATIONARY HAND BRAKE

Push up to unlock and down to lock

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

BOBBY ASTRO 12 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

| | |
|----------------------|-----------------|
| Overall Length | 96.3 cm / 38" |
| Overall Width | 57.5 cm / 22.6" |
| Overall Height | 90.5 cm / 35.6" |
| Folded Length | 69 cm / 27.2" |
| Folded Width | 31 cm / 12.2" |
| Folded Height | 72 cm / 28.3" |
| Seat Depth | 41 cm / 16.1" |
| Seat Width | 48.5 cm / 19.1" |
| Seat Height | 52 cm / 20.5" |
| Arm Rest Height | 69.8 cm / 27.5" |
| Back Rest Height | 37.5 cm / 14.8" |
| Front Wheel Diameter | 15 cm or 6" |
| Back Wheel Diameter | 30.5 cm or 12" |
| Safe Working Load | 115 kg |
| Overall Weight | 12.3 kg |

Check our our full range of detachable wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

