

# BOBBY EVO LIGHTWEIGHT TRAVEL CHAIR



*For illustration purposes only*

**VERMEIREN**

## BOBBY EVO LIGHTWEIGHT TRAVEL CHAIR

### PRODUCT OVERVIEW



### Lightweight Travel Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with detachable foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aircraft grade aluminium lightweight frame
- Designed in Belgium

# VERMEIREN LIGHTWEIGHT TRAVEL WHEELCHAIR

## PRODUCT FEATURES



### FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



### FLIP UP ARM REST

Press button to release lock and flip up arm rest



### DETACHABLE FOOTREST

Pull the lever to release lock and flip up arm rest



### STATIONARY HAND BRAKE

Push up to unlock and down to lock

*For illustration purposes only*

## CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

## BOBBY EVO LIGHTWEIGHT TRAVEL CHAIR

### PRODUCT SPECIFICATIONS

Overall Length	86 cm or 34"
Overall Width (Wheel to Wheel)	53 cm or 21"
Overall Height	83 cm or 33"
Folded Length	66 cm or 26"
Folded Width (Wheel to Wheel)	28 cm or 11"
Folded Height	68.5 cm or 35"
Seat Depth	42 cm or 16.5"
Seat Width	43 cm or 17"
Seat Height (Floor to Seat Base)	48 cm or 19"
Arm Rest Height (Floor to Arm Rest)	68.5 cm or 27"
Back Rest Height	38 cm or 15"
Front Wheel Diameter	15 cm or 6"
Back Wheel Diameter	30.5 cm or 12"
Safe Working Load	120 kg
Overall Weight	15 kg

Check our our full range of lightweight wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

