

FS863-16 LIGHTWEIGHT WHEELCHAIR



FS863-16 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW



Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with flip up arm rest and foot rest for easier patient transfer
- Foldable back rest for easy storage
- Made with beathable woven nylon seat
- Comes with lightweight aluminium frame

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



STATIONARY HAND BRAKE

Push up to unlock and down to lock



ANGLE ADJUSTMENT

Press the button to pull

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

FS863-16 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	106 cm / 41.7"
Overall Width (Wheel to Wheel)	58 cm / 22.8"
Overall Height	89 cm / 35"
Folded Length	76 cm / 29.9"
Folded Width (Wheel to Wheel)	29 cm / 11.4"
Folded Height	64 cm / 25.2"
Seat Depth	39.5 cm / 15.6"
Seat Width	45 cm / 17.7"
Seat Height (Floor to Seat Base)	48 cm / 18.9"
Arm Rest Height (Floor to Arm Rest)	64 cm / 25.2"
Back Rest Height	45 cm / 17.7"
Front Wheel Diameter	15.2 cm / 6"
Back Wheel Diameter	40.5 cm / 16"
Safe Working Load	85 kg
Overall Weight	10.5 kg

Check our our full range of lightweight wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

