

FS870 LIGHTWEIGHT WHEELCHAIR



FS870 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW




Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair for easier patient transfer
- Foldable back rest for easy storage
- Made with beathable woven nylon seat
- Comes with lightweight aluminium frame

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



STATIONARY HAND BRAKE

Push up to unlock and down to lock

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

FS870 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	95 cm or 37.4"
Overall Width (Wheel to Wheel)	58 cm or 22.8"
Overall Height	86 cm or 33.9"
Folded Length	94 cm or 37"
Folded Width (Wheel to Wheel)	29 cm or 11.4"
Folded Height	72.5 cm or 28.5"
Seat Depth	39.5 cm or 15.6"
Seat Width	44.5 cm or 17.5"
Seat Height (Floor to Seat Base)	49 cm or 19.3"
Arm Rest Height (Floor to Arm Rest)	69 cm or 27.2"
Back Rest Height	37 cm or 14.5"
Front Wheel Diameter	17.8 cm or 7"
Back Wheel Diameter	38.1 cm or 15"
Safe Working Load	100 kg
Overall Weight	13 kg

Check our our full range of detachable wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

