

# FT6706 TRAVEL WHEELCHAIR



## FT6706 TRAVEL WHEELCHAIR

### PRODUCT OVERVIEW



### Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with detachable foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aircraft grade aluminium lightweight frame
- Designed in Belgium

*For illustration purposes only*

# LIGHTWEIGHT WHEELCHAIR

## PRODUCT FEATURES



### FLIP UP ARM REST

Flip up arm rest by pushing up



### FLIP UP FOOT REST

Flip up foot rest by pushing up

*For illustration purposes only*

## CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

## FT6706 TRAVEL WHEELCHAIR

### PRODUCT SPECIFICATIONS

Overall Length	83 cm or 32.7"
Overall Width (Wheel to Wheel)	50 cm or 19.6"
Overall Height	87 cm or 34.2"
Folded Length	27 cm or 10.6"
Folded Width (Wheel to Wheel)	48 cm or 18.9"
Folded Height	73 cm or 28.7"
Seat Depth	38 cm or 14.9"
Seat Width	40 cm or 15.7"
Seat Height (Floor to Seat Base)	48 cm or 18.9"
Arm Rest Height (Floor to Arm Rest)	67 cm or 26.3"
Back Rest Height	36 cm or 14.2"
Front Wheel Diameter	15.3 cm or 6"
Back Wheel Diameter	20.3 cm or 8"
Safe Working Load	100 kg
Overall Weight	8 kg

Check our our full range of lightweight wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

