

# JOY RIDER ELECTRICAL WHEELCHAIR



## JOY RIDER ELECTRICAL WHEELCHAIR

### PRODUCT OVERVIEW




### Joy Rider Electrical Wheelchair



Wheelchair sets new standards with its ultra simple, reliable traction drive and intuitive steering system, giving amazing maneuverability. Its unique design creates a revolutionary instant folding frame.

# ELECTRICAL WHEELCHAIR

## PRODUCT FEATURES



### INSERT BATTERY

Always insert one side first by aligning the catch



### SECURE THE BATTERY

Push the lever inward to secure the battery



*For illustration purposes only*

## CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

## JOY RIDER ELECTRICAL WHEELCHAIR

### PRODUCT SPECIFICATIONS

Overall Length	85.5 cm or 33.6"
Overall Width (Wheel to Wheel)	61 cm or 24"
Overall Height	91 cm or 35.8"
Folded Length	111.5 cm or 43.8"
Folded Width	66 cm or 25.9"
Folded Height	32 cm or 12.5"
Seat Depth	33 cm or 12.9"
Seat Width	43.2 cm or 17"
Seat Height	51.8 cm or 20.4"
Arm Rest Height (Floor to Arm Rest)	65 cm or 25.5"
Back Rest Height	44 cm or 17.3"
Front Wheel Diameter	17.7 cm or 7"
Back Wheel Diameter	25.4 cm or 10"
Safe Working Load	115 kg
Overall Weight	23 kg

Check our our full range of electrical wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

