

# KY903 DETACHABLE ELEVATING WHEELCHAIR



*For illustration purposes only*

## KY903 DETACHABLE ELEVATING WHEELCHAIR

### PRODUCT OVERVIEW



### KY903 Detachable Elevating Wheelchair

- Multi-function wheelchair with detachable foot rests and flip up arm rests for easier patient transfer.
- Breathable woven nylon seat
- Lightweight aluminium frame

# ELEVATING WHEELCHAIR

## PRODUCT FEATURES



### FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



### FLIP UP ARM REST

Press button to release lock and flip up arm rest



### DETACHABLE FOOTREST

Pull the lever to release lock and flip the foot rest



### STATIONARY HAND BRAKE

Push up to unlock and down to lock

*For illustration purposes only*

## CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

## KY903 DETACHABLE ELEVATING WHEELCHAIR

### PRODUCT SPECIFICATIONS

|                                     |                  |
|-------------------------------------|------------------|
| Overall Length                      | 129.5 cm or 51"  |
| Overall Width (Wheel to Wheel)      | 62.5 cm or 24.5" |
| Overall Height                      | 94 cm or 37"     |
| Folded Length                       | 86.5 cm or 34"   |
| Folded Width                        | 28 cm or 11"     |
| Folded Height                       | 74 cm or 29"     |
| Seat Depth                          | 46 cm or 18"     |
| Seat Height                         | 50 cm or 19.5"   |
| Arm Rest Height (Floor to Arm Rest) | 70 cm or 27.5"   |
| Back Rest Height                    | 48.5 cm or 19"   |
| Front Wheel Diameter                | 20.5 cm or 8"    |
| Back Wheel Diameter                 | 58.5 cm or 23"   |
| Safe Working Load                   | 100 kg           |
| Overall Weight                      | 16.4 kg          |

Check our our full range of elevating wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

