

ECLK30 KIDS WHEELCHAIR



For illustration purposes only

VERMEIREN

ECLK30 KIDS WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	85 cm / 33.5"
Overall Width (Wheel to Wheel)	50 cm / 19.7"
Overall Height	90 cm / 35.4"
Folded Length	67 cm / 26.3"
Folded Width (Wheel to Wheel)	34 cm / 13.4"
Folded Height	67 cm / 26.3"
Seat Depth	30 cm / 11.8"
Seat Width	31 cm / 12.2"
Seat Height (Floor to Seat Base)	45 cm / 17.7"
Back Rest Height	35 cm / 13.8"
Front Wheel Diameter	12.5 cm / 4.9"
Back Wheel Diameter	50 cm / 19.7"
Safe Working Load	75 kg
Overall Weight	13.6 kg

Check our our full range of detachable wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



ECLK30 KIDS WHEELCHAIR

PRODUCT OVERVIEW



ECLK30 KIDS WHEELCHAIR

- Easy to fold and operate
- Multi-function wheelchair with detachable foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aircraft grade aluminium lightweight frame
- Designed in Belgium

ECLK30 KIDS WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides




ADJUSTABLE ARM REST

Press, hold button and move up or down to adjust or remove



DETACHABLE FOOTREST

Pull the lever to release lock to remove footrest



STATIONARY BRAKES

squeeze and release hand breaks or push up and down wheel breaks to lock and unlock.

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.