

WHEELCHAIRS | USER GUIDE

rainbow care

FS800 TRAVEL WHEELCHAIR



For illustration purposes only

FS800 TRAVEL WHEELCHAIR

PRODUCT OVERVIEW



Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with fold-up foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aluminium lightweight frame
- Pull-Trolley attachment for easy mobility

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FLIP UP ARM REST

Flip up arm rest by pushing up



HANDLE BRAKES

Push up handle lever fully to lock and pull red knob to unlock



FLIP UP FOOT REST

Flip up foot rest by pushing upwards

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates. The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

FS800 TRAVEL WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	85 cm or 33.4"
Overall Width (Wheel to Wheel)	38 cm or 14.9"
Overall Height	88 cm or 34.6"
Folded Length	45 cm or 17.7"
Folded Width (Wheel to Wheel)	33 cm or 12.9"
Folded Height	76 cm or 29.9"
Seat Depth	35 cm or 13.7"
Seat Width	40 cm or 15.7"
Seat Height (Floor to Seat Base)	43 cm or 16.9"
Arm Rest Height (Floor to Arm Rest)	68 cm or 26.7"
Back Rest Height	37 cm or 14.5"
Front Wheel Diameter	15 cm or 5.9"
Back Wheel Diameter	20 cm or 7.8"
Safe Working Load	80 kg
Overall Weight	7 kg

Check our our full range of detachable wheelchair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

