

BATHROOM AIDS | USER GUIDE

rainbow care

FT7210G SHOWER CHAIR



For illustration purposes only

SHOWER CHAIR

FT7210G ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only



Shower Chair

Features a lightweight aluminium frame with adjustable legs to suit every users' height preference, seat handles to grasp on to, and a 12 hole seat with back rest.

SHOWER CHAIR

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

CARE INSTRUCTIONS

- Ensure the shower chair and all fitted parts are secure before use.
- For freestanding models, place the chair on a flat, stable surface and ensure all feet rest firmly on the floor before use.
- For height-adjustable models, ensure all legs are set to the same height and securely locked before sitting.
- Check that rubber feet, where fitted, are secure and are not excessively worn.
- For wall-mounted models, ensure the chair is properly installed on a suitable wall and all mounting fittings are secure. Fully lower the seat into its intended position before use.
- Ensure the seat, backrest, armrests and other fittings, where applicable, are properly positioned and secure before use.
- Do not exceed the safe working load stated for the product.
- When sitting down or standing up, ensure the chair remains stable. Seek assistance if needed.
- Do not stand on the seat or use the chair as a step.
- Clean regularly with a soft cloth and a mild cleaning solution suitable for the product's materials. Follow any material-specific care instructions supplied.
- Rinse away soap residue and dry after use. Avoid abrasive cleaners that may damage the surface.
- Regularly inspect the seat, frame, adjustment points and mounting fittings, where applicable, for looseness, cracks, damage or wear.
- Stop using the chair if any part appears loose, damaged or unstable, and arrange for inspection or replacement before further use.

PRODUCT SPECIFICATIONS

FT7210G

Overall Length	55 cm or 21.7"
Overall Width	50 cm or 19.7"
Overall Height	71 cm - 81 cm 28" - 31.9"
Seat Depth	30.5 cm or 12"
Seat Width (Without Handles)	39 cm or 15.3"
Seat Width (With Handles)	49.5 cm or 19.4"
Seat Height	41 cm - 51 cm 16.1" - 20.0"
Safe Working Load	90 kg
Overall Weight	4 kg

Check out our full range of shower chairs. Review your bathroom environment and consider how tasks, such as transferring someone, are performed before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

