

ECL X1B-24 ECLIPS LIGHTWEIGHT WHEELCHAIR



For illustration purposes only

VERMEIREN

ECL X1B-24 ECLIPS LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW



Lightweight Wheelchair



- Easy to fold and operate
- Multi-function wheelchair for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aircraft grade aluminium lightweight frame
- Designed in Belgium

VERMEIREN LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



STATIONARY HAND BRAKE

Push up to unlock and down to lock

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

ECL X1B-24 ECLIPS LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	111.5 cm or 43.9"
Overall Width (Wheel to Wheel)	70 cm or 27.6"
Overall Height	92 cm or 36.2"
Folded Length	111.5 cm or 43.9"
Folded Width (Wheel to Wheel)	37 cm or 14.6"
Folded Height	71 cm or 28"
Seat Depth	41 cm or 16.1"
Seat Width	42 cm or 16.5" 46 cm or 18.1"
Seat Height (Floor to Seat Base)	47 cm or 18.5"
Arm Rest Height (Floor to Arm Rest)	67 cm or 26.4"
Back Rest Height	40 cm or 15.7"
Front Wheel Diameter	15.2 cm or 6"
Back Wheel Diameter	60.9 cm or 24"
Safe Working Load	90 kg
Overall Weight	14 kg

Check out our full range of wheelchairs. When comparing models, consider the user's environment and how tasks such as transfers will be carried out before deciding

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

