

FS863-12 LIGHTWEIGHT WHEELCHAIR



FS863-12 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW



Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with fixed arm rest and flip up foot rest for easier patient transfer
- Foldable back rest for easy storage
- Made with breathable woven nylon seat
- Comes with lightweight aluminium frame

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES



FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides



STATIONARY HAND BRAKE

Push up to unlock and down to lock



ANGLE ADJUSTMENT

Press the button to pull

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

FS863-12 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	95 cm or 37.4"
Overall Width (Wheel to Wheel)	59 cm or 23.2"
Overall Height	90 cm or 35.4"
Folded Length	79 cm or 31.1"
Folded Width (Wheel to Wheel)	32 cm or 12.6"
Folded Height	66 cm or 26"
Seat Depth	39.5 cm or 15.6"
Seat Width	45 cm or 17.7"
Seat Height (Floor to Seat Base)	48 cm or 18.9"
Arm Rest Height (Floor to Arm Rest)	65 cm or 25.6"
Back Rest Height	45 cm or 17.7"
Front Wheel Diameter	15.2 cm or 6"
Back Wheel Diameter	30.5 cm or 12"
Safe Working Load	80 kg
Overall Weight	10.5 kg

Check out our full range of Lightweight Wheelchairs. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

