WALKING AIDS | USER GUIDE

rainbow care

FS965LH-F6 ROLLATOR



For illustration purposes only

rainbow care

FS965LH-F6 ROLLATOR

PRODUCT OVERVIEW



PRODUCT OVERVIEW



For illustration purposes only

FS965LH-F6 Rollator

The FS965LH-F6 rollator has a flip-up seat, allowing the user to rest when needed. Ensure the brakes are locked before sitting down.

rainbow care

FS965LH-F6 ROLLATOR

PRODUCT FEATURES





HANDLE BRAKES

Squeeze to brake, push down to lock.





ADJUSTABLE HEIGHT

Flip lever up to unlock, adjust height and flip down to lock in place.





FOLDABLE

Pull the center strap upward to fold the rollator.

For illustration purposes only

CARE INSTRUCTIONS

- When using rollator, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each handle is the same before using it.
- Make sure the level of flooring is even before seating down.
- Do not exceed the maximum load.

rainbow care

FS965LH-F6 ROLLATOR

PRODUCT SPECIFICATIONS

Overall Length	64 cm / 25.2”
Overall Width	60.2 cm / 23.7”
Overall Height	81 - 93 cm / 31.9 - 36.6”
Folded Length	64 cm / 25.2”
Folded Width	25 cm / 9.8”
Folded Height	81 cm / 31.9”
Seat Depth	24.5 cm / 9.6”
Seat Width	42.5 cm / 16.7” (Cushion)
Seat Height	52 cm / 20.5” (Floor to Seat)
Front Wheel Diameter	20.5 cm / 8.1”
Back Wheel Diameter	17.5 cm / 6.9”
Safe Working Load	80 kg 5 kg (Basket)
Overall Weight	9 kg

Check out our full range of rollators. When comparing models, consider the user’s environment, mobility needs and required safe working load

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701