WHEELCHAIRS | USER GUIDE

rainbow care

ECL X260 ECLIPS RECLINING WHEELCHAIR



For illustration purposes only

VERMEIREN

rainbow care

ECL X260 ECLIPS RECLINING WHEELCHAIR

PRODUCT OVERVIEW





ECL X260 Eclips Reclining Wheelchair

- Designed for users who require reclining back support, an adjustable headrest and elevating leg rests
- Multi-function wheelchair with detachable foot rests and flip-up armrests for easier patient transfer
- Breathable woven nylon seat
- Lightweight aluminium frame

For illustration purposes only

rainbow care

VERMEIREN RECLINING WHEELCHAIR

PRODUCT FEATURES





RECLINE BACK REST

Press down on the lever to recline the back





FLIP-UP ARM REST

Press button to release lock and flip up arm rest





DETACHABLE FOOT REST

Pull the lever to release lock to remove





ADJUST THE FOOT REST ANGLE

Push the lever to unlock or lock the footrest at the desired angle

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

rainbow care

ECL X260 ECLIPS RECLINING WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	110 cm or 43.3 ”
Overall Width (Wheel to Wheel)	70 cm or 27.6 ”
Overall Height	130 cm / 51.2 ” (With headrest) 104 cm or 40.9 ” (Without headrest)
Seat Depth	48.5 cm / 19.1 ”
Seat Width	45 cm / 17.7 ”
Seat Height	56 cm or 22.0 ”
Reclining Angle	Up to 160°
Arm Rest Height (Floor to Arm Rest)	69 cm or 27.2”
Back Rest Height	67 - 97 cm / 26.4” - 38.2” (With headrest) 54.5 cm or 21.5” (Without headrest)
Front Wheel Diameter	15.2 cm or 6”
Back Wheel Diameter	61 cm or 24”
Safe Working Load	120 kg
Overall Weight	28.5 kg

Check out our full range of Reclining Wheelchairs. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701