WHEELCHAIRS | USER GUIDE

rainbow care

BOBBY 12 LIGHTWEIGHT WHEELCHAIR



For illustration purposes only

rainbow care

BOBBY 12 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW





Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with fold-up foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aluminium lightweight frame

For illustration purposes only

rainbow care

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES





DETACHABLE ARM REST

Pinch and lift the lever to release the armrest lock, then pull upward to detach





CAR BOOT FRIENDLY

Collapse wheelchair into a square to keep in your car boot





FLIP UP FOOT REST

Flip up foot rest by pushing upwards

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

rainbow care

BOBBY 12 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	94 cm or 37”
Overall Width	57.5 cm or 22.6” (Arm rest to Arm rest) 50.5 cm or 19.9” (Wheel to Wheel)
Overall Height	90.5 cm or 35.6”
Folded Length	76 cm or 29.9”
Folded Width (Wheel to Wheel)	33 cm or 13.0”
Folded Height	70.5 cm or 27.8”
Seat Depth	40 cm or 15.7”
Seat Width	47 cm or 18.5”
Seat Height (Floor to Seat Base)	50.5 cm or 19.9”
Arm Rest Height (Floor to Arm Rest)	70 cm or 27.6”
Back Rest Height	40 cm or 15.7”
Front Wheel Diameter	15 cm or 5.9”
Back Wheel Diameter	30.5 cm or 12”
Safe Working Load	115 kg
Overall Weight	12.3 kg

Check out our full range of Lightweight Wheelchairs. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701