WHEELCHAIRS | USER GUIDE

rainbow care

FS800-F5 TRAVEL WHEELCHAIR



For illustration purposes only

rainbow care

FS800-F5 TRAVEL WHEELCHAIR

PRODUCT OVERVIEW





For illustration purposes only

Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with fold-up foot rests for easier patient transfer.
- Foldable back rest for easy storage
- Ventilated and cushioned seat upholstery
- Aluminium lightweight frame
- Pull-Trolley attachment for easy mobility

rainbow care

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES





FLIP UP ARM REST

Flip up arm rest by pushing up





HANDLE BRAKES

Push up handle lever fully to lock and press lower lever to unlock





FLIP UP FOOT REST

Flip up foot rest by pushing upwards

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

rainbow care

FS800-F5 TRAVEL WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	90 cm or 35.4”
Overall Width (Wheel to Wheel)	56.5 cm or 22.2”
Overall Height	93 cm or 36.6”
Folded Length	40 cm or 15.7”
Folded Width (Wheel to Wheel)	47 cm or 18.5”
Folded Height	80.5 cm or 31.7”
Seat Depth	37 cm or 14.6”
Seat Width	48 cm or 18.9”
Seat Height (Floor to Seat Base)	50 cm or 19.7”
Arm Rest Height (Floor to Arm Rest)	74 cm or 29.1”
Back Rest Height	33 cm or 13.0”
Front Wheel Diameter	15 cm or 5.9”
Back Wheel Diameter	28 cm or 11.0”
Safe Working Load	90 kg
Overall Weight	8 kg

Check out our full range of Travel Wheelchair. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701