WHEELCHAIRS | USER GUIDE

rainbow care

# FS863-207-16 LIGHTWEIGHT WHEELCHAIR



*For illustration purposes only*

rainbow care

# FS863-207-16 LIGHTWEIGHT WHEELCHAIR

## PRODUCT OVERVIEW






### Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with flip up arm rest and foot rest for easier patient transfer
- Foldable back rest for easy storage
- Made with breathable woven nylon seat
- Comes with lightweight aluminium frame

*For illustration purposes only*

rainbow care

# LIGHTWEIGHT WHEELCHAIR

## PRODUCT FEATURES





### FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides





### STATIONARY HAND BRAKE

Push up to unlock and down to lock





### FLIP UP ARM REST

Pull the lever to release lock and flip up arm rest

*For illustration purposes only*

## CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

rainbow care

# FS863-207-16 LIGHTWEIGHT WHEELCHAIR

## PRODUCT SPECIFICATIONS

| | |
|---|---|
| Overall Length | 94 cm or 37” |
| Overall Width (Wheel to Wheel) | 57.7 cm or 22.7” |
| Overall Height | 91.4 cm or 36” |
| Folded Length | 78.7 cm or 31” |
| Folded Width (Wheel to Wheel) | 31 cm or 12.2” |
| Folded Height | 64 cm or 25.2” |
| Seat Depth | 40 cm or 15.7” |
| Seat Width | 46 cm or 18.1” |
| Seat Height (Floor to Seat Base) | 48.3 cm or 19” |
| Arm Rest Height (Floor to Arm Rest) | 64 cm or 25.2” |
| Back Rest Height | 44.5 cm or 17.5” |
| Front Wheel Diameter | 15.2 cm or 6” |
| Back Wheel Diameter | 40.6 cm or 16” |
| Safe Working Load | 85 kg |
| Overall Weight | 10.5 kg |

Check out our full range of Lightweight Wheelchairs. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

**Scan the QR code** to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701