WHEELCHAIRS | USER GUIDE

rainbow care

FS863-207-20 LIGHTWEIGHT WHEELCHAIR



For illustration purposes only

rainbow care

FS863-207-20 LIGHTWEIGHT WHEELCHAIR

PRODUCT OVERVIEW





Lightweight Wheelchair

- Easy to fold and operate
- Multi-function wheelchair with flip up arm rest and foot rest for easier patient transfer
- Foldable back rest for easy storage
- Made with breathable woven nylon seat
- Comes with lightweight aluminium frame

For illustration purposes only

rainbow care

LIGHTWEIGHT WHEELCHAIR

PRODUCT FEATURES





FOLDABLE BACKREST

Press on the lever to release and fold back rest downwards for both sides





STATIONARY HAND BRAKE

Push up to unlock and down to lock





FLIP UP ARM REST

Press button to release lock and flip up arm rest

For illustration purposes only

CARE INSTRUCTIONS

- When using the Wheelchair, please pay attention to people and/or items around to avoid collision.
- Lock the wheels before getting in and out of the Wheelchair.
- When getting in and out of the Wheelchair, please do not stand on the footplates.
- The footplates should be folded up / detached beforehand.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.
- When passing through doors or arches, make sure that there is enough room to avoid causing injuries to the user, or the Wheelchair.
- Do not exceed the maximum safe working load.

rainbow care

FS863-207-20 LIGHTWEIGHT WHEELCHAIR

PRODUCT SPECIFICATIONS

Overall Length	104 cm or 40.9”
Overall Width (Wheel to Wheel)	65 cm or 25.6”
Overall Height	90 cm or 35.4”
Folded Length	86 cm or 33.9”
Folded Width (Wheel to Wheel)	37 cm or 14.6”
Folded Height	70 cm or 27.6”
Seat Depth	39.5 cm or 15.6”
Seat Width	46 cm or 18.1”
Seat Height (Floor to Seat Base)	48 cm or 18.9”
Arm Rest Height (Floor to Arm Rest)	67 cm or 26.4”
Back Rest Height	46 cm or 18.1”
Front Wheel Diameter	15.3 cm or 6”
Back Wheel Diameter	50.9 cm or 20”
Safe Working Load	90 kg
Overall Weight	12.5 kg

Check out our full range of Lightweight Wheelchairs. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

6223 0904 • 6931 7425 8501 1878 sales@rainbowcare.com.sg rainbowcare.com.sg
@rainbowcare rainbowcaresg Blk 3007 Ubi Road 1 #02-410 Singapore 408701