

# PA-37WP Wire Alarm

## ALARM MONITOR OPERATION PROCESS



*For illustration purposes only*

### Product Features

- Adjustable sound type and volume
- Included strap
- Water-resistant plastic cover pad
- Battery / Adaptor operated options

### Product Instructions

- The alarm pad is placed on the bed and is connected to an alarm system
- Once the person lay or sits on the pad and applies pressure, the alarm is deactivated
- When the person tries to get off the bed, the pressure is released, activating the alarm and alerting the caregiver or loved one
- Stay alert with our SOS Call Bell

### Product Specifications

Overall Length	51 cm / 20"
Overall Width	76 cm / 30"
Battery	6F22 9 Volts
Operating Frequency	433.92MHz
Working Distance	100 meters (Open Linear Distance)
Volume	Highest: 95-105dB Lowest: 85-95dB
Country Of Origin	Made in China