

KJW322

MOBILE TRANSFER AID



KJW322 MOBILE TRANSFER AID

PRODUCT OVERVIEW



KJW322 Mobile Transfer Aid

- Easy to Use Patient Transfer
- Steel Frame
- Flip Up Foot Rest
- Detachable Bucket
- Safety Belt Included
- 4 Wheels Lock
- Supports up to 110 kg
- Adjustable Height using Cranking System

KJW322 MOBILE TRANSFER AID

PRODUCT FEATURES



SAFETY BUCKLE

Secure the belt for an added layer of safety



RELEASE KNOB

Press the knob to release



HEIGHT ADJUSTMENT LEVER

Precise height adjustment in various applications

CARE INSTRUCTIONS

- Be vigilant of your surroundings to prevent any collisions while using the patient transfer aid.
- Ensure the wheels are locked before getting in and out the patient transfer aid.
- Routinely clean the transfer aid to eliminate dirt, dust, and possible contaminants.
- Follow the weight capacity limits specified by the manufacturer
- Adhere to the specified weight capacity to avoid damage or danger
- Regularly adjust for user comfort—seat cushion and height, etc.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.

KJW322 MOBILE TRANSFER AID

PRODUCT SPECIFICATIONS

Overall Length	75 cm / 29.5"
Overall Width	64 cm / 25.2"
Overall Height	86 - 113 cm / 33.8" - 44.4"
Seat Depth	39.5 cm / 15.6"
Seat Width	53.5 cm / 21.1"
Seat Height	44 - 70 cm / 17.3" - 27.5" (Without Cushion) 49 - 75 cm / 19.2" - 29.5" (With Cushion)
Front Wheel Diameter	10 cm / 3.9"
Back Wheel Diameter	7.5 cm / 3"
Frame Material	Steel
Cushion Cover Material	Foam
Overall Weight	29 kg
Safe Working Load	Up to 110 kg
Country Of Origin	Made in China

Check our our full range of patient hoist. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

