FT3214 WALKING ASSIST





FT3214 ASSIST WALKER AND ROLLATOR

PRODUCT SPECIFICATIONS

Overall Length	72cm / 24.4"
Overall Width	54cm / 21.2"
Overall Height	90cm - 102cm / 35.5"- 40.2" (Adjustable)
Seat Depth	10cm / 3.9"
Seat Width	20cm / 7.8"
Seat Height	25cm / 9.8"
Front Wheel Diameter	14cm / 5.5"
Back Wheel Diameter	15cm / 5.9"

Check our our full range of rollators. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



© 6223 0904 • 6931 7425 © 8501 1878 ≅ sales@rainbowcare.com.sg ⊕ rainbowcare.com.sg ☐ @rainbowcare ☐ rainbowcaresg № Blk 3007 Ubi Road 1 #02-410 Singapore 408701



FT3214 WALKING ASSIST AND ROLLATOR

PRODUCT OVERVIEW



PRODUCT OVERVIEW



For illustration purposes only

FT3214 Walking Assist + Rollator

It's a dual-function mobility aid that supports both walking assistance and rollator use. Its forward-facing design encourages natural posture, while the secure belt system and stable frame provide added safety and confidence.



FT3214 WALKING ASSIST AND ROLLATOR

PRODUCT FEATURES





HANDLE BRAKES FOR BACK WHEELS

Squeeze hand brakes to isolate equipment





FOLDABLE

Pull knob to fold in both Handles for storage

For illustration purposes only

CARE INSTRUCTIONS

- When using rollator, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each handle is the same before using it.
- · Make sure the level of flooring is even before seating down.