

BATHROOM AIDS | USER GUIDE

rainbow care

SHOWER CHAIR



For illustration purposes only

SHOWER CHAIR

FT7500 ALUMINIUM PRODUCT OVERVIEW



FT7500 ALUMINIUM PRODUCT OVERVIEW



Shower Chair

An adjustable shower chair for use in the bathroom. With 4 anti slip legs, it provides a sturdy non slip support even in very wet environments.

SHOWER CHAIR

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull /
push to the level that you
want

For illustration purposes only

PRODUCT SPECIFICATIONS

FT7500

Overall Length 52 cm or 20.5"

Overall Width 52 cm or 20.5"

Overall Height 33 cm - 43 cm
13" - 16.9"

Seat Depth 32 cm or 12.6"

Seat Width 32 cm or 12.6"

Seat Height 40 cm or 15.7"

Safe Working Load 80 kg

Overall Weight 2 kg

CARE INSTRUCTIONS

- When using shower chair, please pay attention to people and/or items around to avoid collision.
- Do not exceed the maximum load.

Check our our full range of shower chair. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

