BATHROOM AIDS | USER GUIDE

rainbow care

STATIONARY COMMODE





| PRODUCT SPECIFICATIONS | FS668B |
|------------------------|------------------|
| Overall Length | 43 cm or 16.9" |
| Overall Width | 37 cm or 14.6" |
| Overall Height | 40 cm or 15.7" |
| Seat Depth | 42 cm or 16.5" |
| Seat Width | 33.5 cm or 13.2" |
| Seat Height | 40 cm or 15.7" |
| Safe Working Load | 70 kg |
| Overall Weight | 2.3 kg |

FS668B PLASTIC PRODUCT OVERVIEW

Stationary Commode

Made of a plastic with functions as a 2 in 1 commode where it is used as a commode chair and a shower chair.



CARE INSTRUCTIONS

- When using stationary commode, please pay attention to people and/or items around to avoid collision.
- Do not exceed the maximum load.

Check our our full range of stationary commodes. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.



© 6223 0904 • 6931 7425
© 8501 1878 ≅ sales@rainbowcare.com.sg ⊕ rainbowcare.com.sg
I @rainbowcare @ rainbowcaresg
© Blk 3007 Ubi Road 1 #02-410 Singapore 408701