

STATIONARY COMMODE



For illustration purposes only

STATIONARY COMMODE

FS894L ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only

FT4100 ALUMINIUM PRODUCT OVERVIEW



For illustration purposes only

STATIONARY COMMODE

FS894L ALUMINIUM



For illustration purposes only

Stationary Commode

Made of lightweight aluminium frame and functions as a 2 in 1 commode where it is used as a commode chair and a shower chair. Foldable for storage.

FT4100 ALUMINIUM



PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want



REMOVE BUCKET

Lift up the seat cover to remove the bucket

For illustration purposes only

| PRODUCT SPECIFICATIONS | FS894L | FT4100 |
|-----------------------------------|--------------------------------|--------------------------------|
| Overall Length | 42 cm or 16.5" | 43 cm or 16.9" |
| Overall Width (Armrest - Armrest) | 55 cm or 21.7" | 46 cm or 18.1" |
| Overall Width (Leg - Leg) | 55 cm or 21.7" | 56 cm or 21.2" |
| Overall Height | 76 cm - 86 cm 29.9" - 33.9" | 60 cm - 70 cm 23.6" - 27.6" |
| Seat Depth | 43 cm or 16.9" | 38 cm or 15" |
| Seat Width (Internal) | 38 cm or 15" | |
| Seat Width (External) | 48 cm or 18.9" | 34 cm or 13.4" |
| Seat Height | 45 cm - 55 cm 17.7" - 21.6" | 43 cm - 53 cm 16.9" - 20.9" |
| Safe Working Load | 80 kg | 100 kg |
| Overall Weight | 5 kg | 4 kg |

CARE INSTRUCTIONS

- When using stationary commode, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each leg is the same before sitting down.
- Do not exceed the maximum load.

Check our our full range of stationary commodes. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

