

# Locomotor Slide Sheet

## TRANSFER AIDS



*For illustration purposes only*

### Product Features

- Patient Transfer Aid
- Reduce Skin Sheering
- Multi Direction Sliding
- Turning a patient for various purposes such as X-rays, surgical procedures, recovery, pressure care, changing dressings, or cleaning
- Sitting a patient up in bed for various reasons such as feeding or providing relief from chest pain
- Moving a patient on or off a bed, trolley, or examination couch without causing painful shearing to the buttocks is crucial for their comfort and well-being

### Product Specifications

Overall Length	122 cm
Overall Width	71 cm / 27.9 "
Material	Nylon
Country Of Origin	Made in United Kingdom