PRODUCT SPECIFICATIONS	FT3502
Overall Length	53.5 cm / 21.1"
Overall Width	46 cm / 18.1"
Overall Height	75 cm - 93 cm 29.5" - 36.6"
Folded Length	11 cm / 4.3"
Folded Width	53 cm / 20.9"
Folded Height	75 cm - 93 cm 29.5" - 36.6"
Safe Working Load	100 kg
Overall Weight	2.4 kg

CARE INSTRUCTIONS

- When using walking frame, please pay attention to people and/or items around to avoid collision.
- · Make sure the level of the height of each leg is the same before using it.
- · Do not exceed the maximum load.

Check our our full range of walking frames. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.



Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

WALKING AIDS | USER GUIDE

rainbow care

FT3502C RECIPROCAL WALKING FRAME





FOLDABLE WALKING FRAME

FT3502A FOLDABLE PRODUCT OVERVIEW



FT3502A FOLDABLE PRODUCT OVERVIEW



Foldable Reciprocal Walking Frame

FT3502A Foldable Reciprocal Walking Frame is suitable for both indoor and outdoor use. The legs are equipped with ferrule tips, providing comfort and security by offering stability and traction.

FOLDABLE WALKING FRAME

PRODUCT FEATURES







Press the button and pull / push to the level that you want





FOLDABLE

Press down the rubber button and fold the side inward

For illustration purposes only

For illustration purposes only