rainbow care

PRODUCT SPECIFICATIONS	KY3320	FT3310
Overall Length	53 cm or 20.9"	50 cm or 19.6"
Overall Width	45 cm or 17.7"	54.4 cm or 21.4"
Overall Height	75 cm - 94 cm 29.5" - 37"	71.5 cm - 89 cm 28.1" - 35"
Folded Length	53 cm or 20.9"	52 cm or 20.5"
Folded Width	10 cm or 3.9"	11 cm or 4.3"
Folded Height	75 cm - 94 cm 29.5" - 37"	71.5 cm - 89 cm 28.1" - 35"
Safe Working Load	100 kg	100 kg
Overall Weight	2.65 kg	2.4 kg

CARE INSTRUCTIONS

- When using walking frame, please pay attention to people and/or items around to avoid collision.
- · Make sure the level of the height of each leg is the same before using it.
- · Do not exceed the maximum load.

Check our our full range of walking frames. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.



Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

WALKING AIDS | USER GUIDE

rainbow care

FOLDABLE RECIPROCAL WALKING FRAME



For illustration purposes only



FOLDABLE RECIPROCAL WALKING FRAME

FT3320 FOLDABLE PRODUCT OVERVIEW



FT3310 FOLDABLE PRODUCT OVERVIEW



FOLDABLE RECIPROCAL WALKING FRAME

FT3320 FOLDABLE





FT3310 FOLDABLE



Foldable Reciprocal Walking Frame

Foldable Reciprocal Walking Frame is a classic design that encompasses several safety features for the overall well-being of its users'.

PRODUCT FEATURES





HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want





FOLDABLE

Press down the rubber button and fold the side inward

For illustration purposes only