

FOLDABLE RECIPROCAL WALKING FRAME



FOLDABLE RECIPROCAL WALKING FRAME

FT3320 FOLDABLE PRODUCT OVERVIEW



For illustration purposes only

FT3310 FOLDABLE PRODUCT OVERVIEW



For illustration purposes only

FOLDABLE RECIPROCAL WALKING FRAME

FT3320 FOLDABLE



FT3310 FOLDABLE



For illustration purposes only

Foldable Reciprocal Walking Frame

Foldable Reciprocal Walking Frame is a classic design that encompasses several safety features for the overall well-being of its users'.

PRODUCT FEATURES



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want



FOLDABLE

Press down the rubber button and fold the side inward

For illustration purposes only

| PRODUCT SPECIFICATIONS | FT3320 | FT3310 |
|------------------------|--------------------------------|---------------------------------|
| Overall Length | 48 cm / 18.9 " | 47 cm / 18.5 " |
| Overall Width | 57 cm / 22.4 " | 57 cm / 22.4 " |
| Overall Height | 77 cm / 30.3 " 94 cm / 37 " | 72 cm / 28.3 " 89 cm / 35 " |
| Folded Length | 54 cm / 21.3 " | 10 cm / 3.9 " |
| Folded Width | 10 cm / 3.9 " | 60 cm / 23.6 " |
| Folded Height | 75 cm - 94 cm 29.5" - 37" | 71.5 cm / 28.1" 89 cm / 35 " |
| Safe Working Load | 100 kg | 100 kg |
| Overall Weight | 2.65 kg | 2.4 kg |

CARE INSTRUCTIONS

- When using walking frame, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of each leg is the same before using it.
- Do not exceed the maximum load.

Check our our full range of walking frames. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

