

FS9112L WALKING STICK WITH SEAT



For illustration purposes only


PRODUCT SPECIFICATIONS

FS9112L

| | |
|------------------------|---------------------------------------------|
| Overall Height | 81 cm - 89 cm or 31.9" - 35" |
| Base (Length / Width) | 30 cm or 11.8" / 21 cm or 8.3" |
| Seat (Height / Width) | 22 cm or 8.7" / 46 - 52 cm or 18.1" - 20.5" |
| Safe Working Load | 70 kg |
| Overall Weight | 1 kg |

PRODUCT FEATURES

It is a great mobility aids for individuals that don't require the full assistance that a walker or rollator can provide. In most cases, a cane can assist the user in getting up and supporting weight on one or both of their legs.



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

For illustration purposes only

CARE INSTRUCTIONS

- When using walking stick, please pay attention to people and/or items around to avoid collision.
- Make sure the level of the height of leg is the same before using it.
- Make sure notice the wear and tear of rubber stopper and replace when it is needed.
- Do not exceed the maximum load.

Check our our full range of walking sticks. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

☎ 6223 0904 • 6931 7425 ☎ 8501 1878 ✉ sales@rainbowcare.com.sg 🌐 rainbowcare.com.sg

📱 @rainbowcare 📷 rainbowcaresg 📍 Blk 3007 Ubi Road 1 #02-410 Singapore 408701

