

FT2644 WALKING STICK WITH SEAT



For illustration purposes only



PRODUCT SPECIFICATIONS

FT2644

Overall Height	87 cm - 97 cm or 34.3" - 38.2"
Seat Depth / Width	20.5 cm or 8.1" / 24.2 cm or 9.5"
Seat Height	49 cm - 59 cm or 19.3" - 23.2"
Safe Working Load	80 kg
Overall Weight	0.4 kg

PRODUCT FEATURES

It is a great mobility aids for individuals that don't require the full assistance that a walker or rollator can provide. In most cases, a cane can assist the user in getting up and supporting weight on one or both of their legs.



HEIGHT ADJUSTMENT

Press the button and pull / push to the level that you want

For illustration purposes only

CARE INSTRUCTIONS

- Make sure the level of the height of leg is the same before using it.
- Make sure notice the wear and tear of rubber stopper and replace when it is needed.
- Do not exceed the maximum load.

Check our our full range of walking sticks. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

