

# ONTARIO 2 GERIATRIC CHAIR



*For illustration purposes only*

## ONTARIO 2 GERIATRIC CHAIR

### PRODUCT FEATURES



#### ELECTRONICALLY CONTROLLED

Use the intuitive electronic remote to adjust the chair with ease.



#### SMOOTH RECLINE

Fully recline the chair for maximum comfort by using the remote.



#### POWER LIFT ASSIST

Activate the power lift mechanism to help stand up or sit down safely.

*For illustration purposes only*

### CARE INSTRUCTIONS

- Keep the chair on a level surface at all times to ensure stability
- Do not exceed the maximum load specified by the manufacturer
- Regularly check/tighten all screws and moving parts to maintain proper function
- Use the safety belt to prevent the user from sliding out of the chair
- Keep the upholstery clean by wiping with a damp cloth; avoid harsh chemicals
- Ensure that the chair is positioned with sufficient space around it to prevent damage
- Always use the lift and recline functions gently to avoid wear and tear
- Check the power cord and battery regularly for damage or loose connections.

## ONTARIO 2 GERIATRIC CHAIR

### PRODUCT OVERVIEW



*For illustration purposes only*

Normandie Geriatric Chair is a Modern Take on a Classic, designed for the user's needs and comfort. Perfect in its contours, it features a maximum backrest recline of 130° with accompanying extended Foot Rest, Detachable Table, comfortable high-density foam upholstery and generous seating space.

It comes with Seat Belt and Wheel Locks on all four castors to ensure it stays in a secure position. Arm Rests are height adjustable to suit unique individuals' preference. Head Rest can also be easily removed and adjusted accordingly.

## PRODUCT SPECIFICATIONS

### ONTARIO 2 GERIATRIC CHAIR

Overall Length	96 cm / 37.8"
Overall Width	84.5 cm / 33.3"
Overall Height	114.5 cm / 45.1"
Seat Depth	53 cm / 20.9"
Seat Width	48 cm / 18.9"
Seat Height	50 cm / 19.7"
Reclining Position	+9° / +9° / -18°
Frame Material	Steel
Seat Material	PVC Leather
Overall Weight	53 kg
Safe Working Load	150 kg
Designed In	Belgium
Country Of Origin	Made In India

Check out our full range of Reclining Geriatric Chairs. When making comparisons, our best suggestion would be for you to think about your environment and consider how certain task like transferring someone is done before deciding.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

