

KJW325

MOBILE TRANSFER AID



KJW325 MOBILE TRANSFER AID

PRODUCT OVERVIEW



KJW325 Mobile Transfer Aid

- Easy to Use Patient Transfer
- Steel Frame
- Flip Up Foot Rest
- Detachable Bucket
- Safety Belt Included
- 4 Wheels Lock
- Supports up to 100 kg
- Adjustable Height

KJW325 MOBILE TRANSFER AID

PRODUCT FEATURES



SAFETY BUCKLE

Secure the belt for an added layer of safety



RELEASE KNOB

Pull the knob to release



LOCKABLE WHEELS

Step on the ON to lock the wheels and OFF to unlock

CARE INSTRUCTIONS

- Be vigilant of your surroundings to prevent any collisions while using the patient transfer aid.
- Ensure the wheels are locked before getting in and out the patient transfer aid.
- Routinely clean the transfer aid to eliminate dirt, dust, and possible contaminants.
- Follow the weight capacity limits specified by the manufacturer
- Adhere to the specified weight capacity to avoid damage or danger
- Regularly adjust for user comfort—seat cushion and height, etc.
- Pay attention to the shift of the Center of Gravity when on elevated grounds.

KJW325 MOBILE TRANSFER AID

PRODUCT SPECIFICATIONS

Overall Length	70 cm / 27.6"
Overall Width	54 cm / 21.3"
Overall Height	88.5 - 106 cm / 34.8" - 41.7"
Seat Depth	39.5 cm / 15.6"
Seat Width	48 cm / 18.9"
Seat Height	45 - 62 cm / 17.7" - 24.4" (Without Cushion) 50 - 67 cm / 19.7" - 26.4" (With Cushion)
Front Wheel Diameter	13 cm / 5.1"
Back Wheel Diameter	7.5 cm / 3"
Frame Material	Steel
Cushion Cover Material	Foam
Overall Weight	18 kg
Safe Working Load	Up to 100 kg
Country Of Origin	Made in China

Check out our full range of mobile transfer aid. When making comparisons, our best suggestion would be to review your home environment and consider how certain tasks, such as transferring someone, are performed.

Scan the QR code to find out more about the product. To scan, use your camera app or a QR code reader on your device.

