

Wired Bed Alarm Pad

WITH PRESSURE SENSOR



For illustration purposes only

Product Features

- Pressure Sensor Alarm
- Water Resistant Cover
- Suitable for dementia patients or patients who are at risk of falling
- Adjustable volume
- 2 different sounds to choose from

Product Instructions

- The alarm pad is placed on the bed and is connected to an alarm system
- The alarm may be deactivated manually by pressing the button on the device, or once the patient lays or sits back down
- Suitable for dementia patients or patients who are at risk of falling
- When the person tries to get off the bed, the pressure is released, activating the alarm and alerting the caregiver promptly

Product Specifications

Overall Length	51 cm / 20"
Overall Width	76 cm / 30"
Country of Origin	Made in China
Battery	6F22 9 Volts
Country Of Origin	Made in China